May 2025

*Kwik*Blot Cleaning Solution

(Cat.# R2008)

Usage

This buffer is formulated to clean ***Kwik*Blot** automatic blot processor after each use. It prevents rust and mold accumulation on internal metal valves and tubings. This is a crucial part of device maintenance, and we provide a free warranty service if the user uses this solution.

Component

100x ***Kwik*Blot** Cleaning Solution: 40mL

Procedure

1. Dilute the 100x solution in water to make a 1x working solution.
2. Add 40mL the 1x solution on the trays (excluding #1 and #2 trays).
3. Follow the on-screen command of ***Kwik*Blot** phone app to clean the device.